

# ZIMMER® TISSUE BANK DERMATOME



## CONTROL THE POWER

With more than 50 years of experience in skin graft instrument design and manufacturing, Zimmer is pleased to introduce the *Zimmer* Tissue Bank Dermatome. This dermatome is designed to meet the special needs of tissue banks in retrieving cadaver skin. It features a wide range of cutting thickness options from 0 to .060in. (0 to 1.52 mm) with a cutting speed of 5,000 to 6,000 cycles per minute. A heavy-duty motor, along with special stainless steel cutting blades, gives the *Zimmer* Tissue Bank Dermatome the ability to retrieve quality tissue with minimal effort.

## BETTER CONTROL, LESS FATIGUE

Zimmer has incorporated extra power and durability in the motor for the *Zimmer* Tissue Bank Dermatome without the extra weight that can compromise maneuverability. The lightweight unit is durable yet easy to hold and control with less fatigue for the operator. For easy maintenance, the motor never requires lubrication.

## STEAM AUTOCLAVABLE = LESS DOWNTIME

The *Zimmer* Tissue Bank Dermatome can be steam autoclaved for minimal downtime between procedures. An optional autoclave case conveniently holds all components for easy storage and transportation.

## QUALITY CUT FROM PRECISION BLADES

Honed from stainless steel to maximum sharpness, these blades require fewer changes during procedures. Pre-sterilized for added convenience. Blades provide consistent quality grafts, even cutting and cost efficiency.

## ONE-YEAR WARRANTY PROTECTION

All parts and labor are guaranteed for one year against normal wear and use. (Annual calibration check recommended.) Servicing normally requires only two to three days. Loaner units are available to minimize any disruption in tissue bank schedules or procedures.

## CONTROL SWITCH LOCATED FOR CONVENIENCE

The ON/OFF Switch for the Tissue Bank is conveniently located on the top of the Dermatome handle for easy activation. A SAFE position helps prevent activation when the Dermatome is not in use, or during blade changes.



**zimmer**  
Confidence in your hands™

# THE QUALITY YOU EXPECT...THE POWER YOU NEED

## ORDER INFORMATION

Prod. No.	Description
8831-000-00	Tissue Bank Dermatome with Screwdriver and Width Plates (3 in. and 4 in.)
8831-003-00	Sterilizing Case
8831-004-00	3 in. (8cm) Width Plate
8831-005-00	4 in. (10cm) Width Plate
8831-000-10	Blades, Box of 10, Sterile
8821-006-00	Power Supply

## ZIMMER TISSUE BANK DERMATOME

- Wide range of graft thickness options, from 0 to .060 inches.
- Lightweight and streamlined for easy maneuverability and less hand fatigue.
- Can be steam autoclaved for convenience and fast turnaround.
- Blades can be changed quickly with less chance of operator injury. No lubrication required.
- Width plates allow option of grafts in either 3 in. or 4 in. widths.
- Power Supply offers options for 110, 120 or 240 volts for any power source.
- 15-foot power cord keeps power supply away from sterile field.

